

Bode SecurSwab DNA Collection Systems

ENHANCED DNA RECOVERY



Bode Cellmark Forensics' **SecurSwab DNA Collection Systems** integrate the collection, labeling, protection, preservation, transportation, processing and storage of DNA evidence.

Bode SecurSwab DNA Collection Systems

Bode SecurSwab DNA collection devices are complete collection systems designed to improve DNA recovery from evidentiary samples.



More DNA

In Bode comparisons, Bode SecurSwab collection devices were shown to recover up to 30% more DNA than other products.¹



Easy to Use

The unique SecurSwab design allows for collection and storage in one single device, eliminating the need to purchase and assemble a swab box for sample storage.



Saves Time

There is no need for time consuming air drying of the swab that may increase the risk of contamination. The integrated desiccants begin drying a collected sample in the protective tube immediately at the crime scene.

Number One-ranked Swab

The Bode SecurSwab DUO-V swab system was the number one-ranked swab for both small and large volume DNA recovery, in the Swab Collection Study conducted by the National Forensic Science Technology Center and funded through the National Institute of Justice.²



DNA degradation begins the moment a sample is collected. **Keep your valuable evidence secure with the Bode BioSafe Swab.**



New Product Spotlight: **Bode BioSafe Swab**

- Point-of-collection preservation for DNA evidence
- Improved DNA Yield and Stability
 - The Bode BioSafe Swab recovered significantly more DNA than a cotton swab stored in a swab box:
 - 76x more DNA from blood stains on a wood surface*
 - 12x more DNA from blood stains on a shoe*
 - 50x more DNA from saliva stains*

*Comparisons were based on the average amount of DNA recovered from the small target in a controlled internal accelerated stability study.

Bode SecurSwab Products



Bode SecurSwab 2 DNA Collection System

The Bode SecurSwab 2 is a complete collection system designed with a protective tube, barcode labels, and integrated desiccants.

This design provides a quick-drying system to protect and preserve evidence, simplifies chain-of-custody, ensures the shaft won't break during collection, and minimizes contamination.



Bode SecurSwab DUO-V DNA Collector

The Bode SecurSwab DUO-V is a complete collection system with two cotton-tipped swabs for wet/dry collection and two integrated desiccants for faster drying. A fully-saturated swab dries in less than six hours.



Bode SecurSwab 2 DNA Collection Kit

Bode SecurSwab crime scene collection kits are fully customizable to contain exactly what an agency or investigator needs to collect and document a DNA sample.



Bode IFP (Intelligence Field Pack) DNA Collection System

Conveniently packed as 6 collectors per pouch, the Bode IFP is designed for fast and efficient collection while maintaining chain-of-custody.

The Bode Cellmark Advantage for DNA Collection Systems

Trust	Bode Cellmark has delivered more than five million collection devices worldwide.
Experience	Decades of experience make Bode Cellmark one of the world's most trusted companies for supplying collection devices for the forensic DNA market.
Innovation	Bode Cellmark continues to develop products and services that enhance our offerings and make us a full-service solutions provider.
Technical Expertise	Bode Cellmark is at the forefront of research and development of DNA testing and provides a high level of quality and service.
Customer Service	Bode Cellmark uses a dedicated technical support team to provide unparalleled customer service for the support and consultation of clients.

Product	Quantity	Catalog Number
Bode SecurSwab 2	Case of 100 collectors	P13D04
Bode SecurSwab 2 with barcodes	Case of 100 collectors	P13D05
Bode SecurSwab 2 collection kit	Case of 25 kits	P13D77
Bode SecurSwab DUO-V collectors with barcodes	Case of 100 collectors	P10D15
Bode IFP collectors	Case of 60 collectors Case of 240 collectors	P10D10 P10D14
Bode DNA collector transport pouch	Case of 100 pouches	P10D43
3 mL water vials	Case of 100 vials	P08D75
Bode BioSafe Swab	Case of 100 collectors	P18D01
Bode BioSafe Swab with barcodes	Case of 100 collectors	P18D02

Bode SecurSwab DNA Collection System Ordering Chart

Answers with confidence and accuracy

Request a sample or speak to a Bode Cellmark representative about DNA Collection offerings:

866-263-3443 or **703-646-9740** or **BTBodeService@LabCorp.com**



References

1. Internal data. Bode Technology, Lorton, Va, Drying Time and its Effect on DNA Yield.
2. National Forensic Science Technology Center. Swab Collection Study. http://www.nfstc.org/?dl_id=267. Published April 24, 2012. Accessed July 22, 2015.
3. Internal data. Bode Cellmark Forensics, Lorton, VA. BioSafe Swab Mock Crime Scene-Modified Condensed Report.

10430 Furnace Road, Suite 107
Lorton, VA 22079